

Guard Column

Instruction Manual

1. Precautions for Column

- This guard cartridge column is packed with packing materials of various functional groups.
- Please use this guard column with the same packing material as analytical column.
- The correct direction of the solvent flow is indicated by an arrow on the column identification label.
- Please read instruction manual of the analytical column before use.
- Early replacement of the guard column is recommended to prevent contamination of analytical column.
- Please be aware that specifications are subject to change without prior notification.

Columns are shipped under highly controlled conditions. However, if you should find any defect, please contact your dealer or Shinwa. Note that Shinwa does not warrant the product against column life or deterioration caused by the failure to follow the above instructions.